

910-0019 Cleaning Solution 1**50 ml**

Contents	50 ml Cleaning Solution 1 The solution contains sodium hypochlorite intended for cleaning of flow-through sample compartment (pH 11.5-13.5). EINECS no. for sodium hypochlorite is 231-668-3
Application	Cleaning Solution 1 is intended for cleaning of the XyctoMatic® internal flow system and flow cell.
Principle	Cleaning Solution 1 contains active oxidation agents, which can dissolve lipids and proteins that might adhere to the windows of the sample compartment
Use	Please refer to the XyctoMatic® User Guide and/or Cleaning Application Note.
Storage	Store at room temperature and protect from light.
Stability	For sealed containers refer to the expiry date on the label of the container. After breaking the seal of the container Cleaning Solution 1 expires after 12 months (room temperature). During this period the cap must be mounted on the container while not in use.
Safety information	Cleaning Solution 1 contains active components, which can be irritating to skin and eyes. Please refer to the Safety Data Sheet (SDS), which is available for (Download) online.
Disposal	After use, dispose of the bottle following national or local laws and regulations, according to the nature of its contents.